



HDX Resuscitation cases

The HDX are surely cases designed for special uses, because they are built to resist impact, corrosive substances, water, dust and weather conditions.

The internal part can be set up with diverse solutions:

- internal padding in previously drilled soft polypropylene with a matrix of small cubes that have a dimension of 1cm², so to obtain the same precise shape of the object to be transported and contained in the case
- internally arranged for inserting metallic supports to allow setting up panels or masks for the fixation of devices or instruments

The HDX cases undergo severe testing and have obtained the following certifications:

IP67 CERTIFICATE
DEF-STAN 81-41 part 3/ issue 4
STANAG 4280 Ed. 2

Dimensions: 520x435x230 mm

HDX empty
HDX complete
HDX complete with defibrillator

Code CB06042C
Code CB06043C
Code CB06044C





**Lienz
Resuscitation kit**

In order to achieve top quality it is necessary to combine a man's love for his job and for detail. It is through this artisanal performance that we can reach the perfection Spencer's cases stand out for. Heavy testing ensures that each and every case responds to high standards of quality. At the end of the day the Spencer brand is a worldwide symbol for high quality. Cases conceived and developed by critical care doctors to optimise location of products and time management during emergencies.

Dimensions Lienz: 520x380x220 mm
Weight Lienz complete: 6,2 kg

Lienz empty
Lienz complete adult
Lienz complete paediatric

Code CB06002A
Code CB06018A
Code CB06019A





Max 911
Resuscitation case

Made of strong tripolymer plastics which are shiny smooth and easy to clean. A neoprene seal is fitted to make them dustproof. A wall bracket secure them on any house wall or booth. The closing security is given by two swivel catches.

Includes:

- 1 Oxygen tank 0,46 lt with pressure reducer
- 1 Oxygen mask low flow adult
- 1 Co B-life PVC manual resuscitator adult without overpressure valve
- 1 Cuffed Mask PVC mask size 3
- 1 Cuffed Mask PVC mask size 4
- 1 Cuffed Mask PVC mask size 5
- 1 Moulder jaw spreader
- 1 Guedel airway non sterile size 3
- 1 Guedel airway non sterile size 4
- 1 Guedel airway non sterile size 5

Dimensions: 400x300x150 mm
Weight complete: 5,8 kg

Max 911 complete

Code MX09300G





Spencer 2
Resuscitation case

Made of strong tripolymer plastics which are shiny smooth and easy to clean. A neoprene seal is fitted to make them dustproof. A wall bracket secure them on any house wall or booth. The closing security is given by two swivel catches.

Includes:

- 1 Oxygen tank 0,46 lt with pressure reducer
- 1 Anti-crushing connective tubing
- 1 Pro B-life PVC manual resuscitator adult
- 1 Cuffed Mask PVC mask size 3
- 1 Moulder jaw spreader
- 1 Collin stainless steel tongue forceps
- 1 Collin techno polymer tongue forceps
- 1 Guedel airway non sterile size 1
- 1 Guedel airway non sterile size 2
- 1 Guedel airway non sterile size 3
- 1 Guedel airway non sterile size 4

Dimensions: 385x250x135 mm
Weight complete: 4,750 kg

Spencer 2 empty
Spencer 2 complete

Code MX09010G
Code RM80000B





Matrix 100
Resuscitation case

Made of semi-rigid thermo-moulded material. Inside, the nylon compartments let you create ideal spaces for other small life-saving objects. Thanks to its strong and waterproof material and to its zipper closure, the inside objects are protected from bumps and isolated from high temperatures.

Includes:

- 1 Co B-life PVC manual resuscitator adult without overpressure valve
- Guedel airway non sterile size 3
- Guedel airway non sterile size 4
- Guedel airway non sterile size 5
- 1 Moulder jaw spreader
- 1 Young tongue forceps
- 1 Cuffed Mask PVC mask size 4

Dimensions: 425x331x128 mm

Matrix 100 empty
Matrix 100 complete

Code BS05010A
Code RM90000E





REM
Resuscitation cases

With the REM collection Spencer presents a new generation of lightweight cases, conceived, designed, shaped and realised using the latest techniques. The casing is made in aluminium 100%, a material that has changed the world and with which we can still make a lot of dream ideas come true. A line thought for those who love dynamic as well as practical accessories.

Includes:

- 1 Co B-life PVC manual resuscitator adult without overpressure valve
- Guedel airway non sterile size 3
- Guedel airway non sterile size 4
- Guedel airway non sterile size 5
- 1 Moulder jaw spreader
- 1 Young tongue forceps
- 1 Cuffed Mask PVC mask size 4

Dimensions: 455x335x130 mm

REM Case empty
REM Case complete
REM Case complete with defibrillator

Code CB06090A
Code CB06092A
Code CB06091A

