



R-OX
Oxygen tank holder

The new thing about R-OX is its extraordinary design, its material and its easy use. The structure, made of welded tubes of aluminium sheets, is very light and at the same time it guarantees a shiny finishing very easy to clean. The position of the handles are unique and assure a comfortable handgrip even for left-handers.

The rear support allows to keep the trolley and tank in an inclined position. All the weight is completely loaded upon the wheels, the handle is only used to drive the trolley avoiding the operator to carry all the weight of the trolley and tank. The straps that hold the tank also have rubber bands and adjustable polypropylene belts to guarantee a safe attachment, protection and stability during transport. Two large wheels with sphere cushions guarantee easy handling and no noise on rough terrain.

Thanks to a generous loading platform and to the innovative supports, it can contain different sizes of tanks. The supports that have nitrylic integrated elements can hold 3, 5 or 10 lt tanks. When not in use, it can be easily folded and stored.

- | | |
|---|---------------|
| R-OX oxygen tank holder/trolley | Code OX10400E |
| R-OX with oxygen tank 5 lt, RAP flowmeter with humidifier and pressure reducer | Code OX10401E |
| R-OX with oxygen tank 10 lt, RAP flowmeter with humidifier and pressure reducer | Code OX10402E |





REM Trolley
Resuscitation trolley

In an era of complete medical device outfits and of large mobility, there is a need to produce cases that can be reliable and efficient companions for you. This is where the REM trolley springs from. Very spacious and accurately designed in the details, the REM is internally structured with two compartments, one of which is a pre-modelled mould in ABS specifically shaped to contain and solidise the most sophisticated devices.

The two shells are integrated with atraumatic divisors and flaps. The body of the case is made with a cold molding process from a 12/10 sheet of anodised aluminium, the corners are protected with an application of tailored alloy sheet, and metallic reinforcements on the profile seal will guarantee perfect tightness.

REM is equipped with a retractile handle and two extra fixed handles accurately positioned; the rubber wheels have a wide diameter and are particularly efficient and robust.

Dimensions: 495x335x220 mm

REM Trolley empty
REM Trolley complete

Code CB06100A
Code CB06101A

