

SKID EVACUATION CHAIR

Transport of disabled persons in evacuation situations

Leader in the market for quality, innovation and image, Spencer represents the most qualified technological reference in terms of evacuation.

In fact, over the years, Spencer evacuation chairs have become synonymous of high quality and prestige all over the world of emergency: a professional benchmark able to combine excellent tradition with technology.

The wide range of models right now established in the market are proof of their indisputable success.

Skid Series

The most relevant technical/aesthetic novelty of this article is the adoption of passive caterpillar belts which effortlessly glide downstairs, with an elevated level of security.

On landing, the constructive geometry of the unit transfers all weight on the main wheels to allow turning and movement.


Moreover, the Pro Skid version is equipped with front AND rear telescopic handles which allow upstairs transport.

A careful examination of the whole chair will reveal a seat that has been very carefully studied; aesthetically and in its conceptually free styling.

The aggressivity has been achieved by the mixing of technical requirements. For example the tortuous conformation of the welded frame or the obvious originality of the seating posture.

The precocity of the aesthetics is underlined by the distribution of the angles, the backrest that can be telescopically inclined, the adjustable headrest and the handles for transportation that can also be detached from the chair itself.

Innovative design and the maximum of security possible when in use, make the new Skid Series an absolute novelty within the environment of evacuations and transport.



www.evacuationchair.net

SPENCER[®]

Spencer Italia S.r.l.
Strada Cavi, 7
43044 Collecchio (Parma) Italy
Tel. +39 0521 541111
Fax. +39 0521 541222
Email: export@spencer.it
www.spencer.it